



THIS IS RECOVERY

THE FULLY INTEGRATED, ACTIVE COLD AND COMPRESSION SYSTEM

The ultimate in integrated cold and compression therapy and clinically proven outcomes. The powerful Game Ready® System and complete range of anatomically designed, dual-action wraps allow for optimal therapeutic benefits.



Temperature, pressure, and time status indicators provide clear information, making it easy to monitor and manage synergistic compression and cold treatments with simple buttons and dials



Adjustable set-point temperature control provides consistent temperature across the wrap and throughout treatment, helping to improve patient comfort



Programs are designed to promote patient compliance and reduce risk of injury, such as frostbite, by automating alternating treatment and rest times




G A M E  R E A D Y®

GAME READY COLD WORKS FASTER, PENETRATES DEEPER, AND LASTS LONGER

Game Ready's active compression squeezes greater benefits out of cold therapy alone. Circumferential wraps envelop the injury or surgery site and pneumatic compression works to conform the wrap to the contours of the body, allowing for better surface contact and more effective cooling.

Head to head infrared thermal imagery shows that Game Ready provides faster, deeper, and longer-lasting cold than other cryotherapy products.



Thermal Imagery: Side view shown. First skin temperature measurement immediately after completing 30 minutes of treatment, subsequent images at 10 and 30 minutes post-treatment. All devices applied per manufacturer's instructions and all measurements taken from uniform distance. All devices applied on different days allowing full recovery for test subject.

PRESCRIBED AND USED BY THOUSANDS OF CLINICIANS, APPRECIATED BY THOUSANDS OF SATISFIED PATIENTS

Game Ready can be found in the most prestigious orthopedic centers, physical therapy clinics, and athletic training facilities. In fact, the treatment system is relied on by Olympic athletes, special military lig forces and countless teams from virtually every professional sport around the world.

HELP PATIENTS REDUCE PAIN & OPIOIDS

Compared to traditional RICE application (REST-ICE-COMPRESSION-ELEVATION), Game Ready delivers greater therapeutic benefits:

- Reduces patient opioid consumption and may allow for earlier discontinuation of opioid medications in certain patients ^{1,2,3}
- Increased post-operative range of motion in ACL patients ³
- May increase patient satisfaction with the recovery process ^{1,4}
- Decreases pain ²

ATX DUAL-ACTION WRAPS

Game Ready ATX® (ACTIVE TEMPERATURE EXCHANGE) wraps feature two internal chambers—one for air, the other for ice water. The patented construction delivers a no-mess integration of adjustable compression and cold therapies. Unlike one-design-fits-all products, Game Ready's comprehensive line of anatomically designed wraps allows for optimal comfort and therapeutic benefits.



DISCOVER THE GAME READY DIFFERENCE

To learn more about the full-service Game Ready Rental Program, visit us at gameready.com or contact Patient Care at **1.888.426.3732**

GAME  READY®

HUMAN INSPIRED. TECHNOLOGY DRIVEN.™

1. Su EP, Perna M, Boettner F, Mayman DJ, Gerlinger T, Barsoum W, Randolph J, Lee G. A prospective, multi-center randomized trial to evaluate the efficacy of a cryopneumatic device on total knee arthroplasty recovery. J Bone & Joint Surg. November 2012; 94-B; SuppA:153-6. 2. Waterman W, Walker JJ, Swains C, Shortt M, Todd MD, Machen SM, Owens BD. (2012). The efficacy of combined cryotherapy and compression compared with cryotherapy alone following anterior cruciate ligament reconstruction. J Knee Surg. 2012;25(02): 155-160. 3. Murgler J, Cassard X. Cryotherapy with dynamic intermittent compression for analgesia after anterior cruciate ligament reconstruction. Preliminary study. Orthop Traumatol Surg Res. 2014;100(3):309-312. 4. Leegwater NC, Jore HW, Brohet R, Nolte P. Cryocompression therapy after elective arthroplasty of the hip. Hip Int. 2012;22(05):527-533.